

B-LITE®
PRODUCT CATALOGUE
N° 5
1 JUNE 2019

B · L I T E

Made in
Germany

The New Power-Brand of
POLYTECH

IMPLANT SIZE SELECTION FOR DIMENSIONAL PLANNING



ISO 14607:2018-04

B-Lite® implants are available in round or anatomical shapes and in a variety of projections. The large variety of sizes combined with the systematic approach to surgical planning ensure that the correct implant is available for every patient.

Base measurements A and B provide the width and height of the implant base respectively. All implants in the catalog have a round base and therefore A and B measurements are equivalent.

Projection measurement C provides the maximum projection of the implant.







Lower coverage measurement D provides the distance between the highest point of projection and caudal implant edge. It creates a correlation between the correct placement of the Nipple Areolar Complex (NAC) and the inframammary fold (IMF).



Note: B and C measurements have changed positions versus prior editions of this catalogue and the B-Lite® box.

ROUND RR

textured

PROFILE	BASE	PROJECTION
MODERATE		
HIGH		
EXTRA HIGH		







ROUND - HIGH PROJECTION					
ORDER N°	DIMENSIONS (mm)				VOLUME (ml)
	A	B	C	D	
20626-165LW25	93	93	39	64	165
20626-190LW25	98	98	41	67	190
20626-215LW25	103	103	43	71	215
20626-250LW25	108	108	44	74	250
20626-285LW25	113	113	46	77	285
20626-325LW25	118	118	48	80	325
20626-365LW25	123	123	49	83	365
20626-410LW25	128	128	51	87	410
20626-460LW25	133	133	53	90	460
20626-510LW25	138	138	54	93	510
20626-565LW25	143	143	56	96	565
20626-625LW25	148	148	58	100	625
20626-690LW25	153	153	59	103	690
20626-760LW25	158	158	61	106	760

ROUND - MODERATE PROJECTION					
ORDER N°	DIMENSIONS (mm)				VOLUME (ml)
	A	B	C	D	
20625-130LW25	93	93	32	59	130
20625-150LW25	98	98	34	62	150
20625-175LW25	103	103	35	65	175
20625-200LW25	108	108	36	68	200
20625-230LW25	113	113	37	71	230
20625-260LW25	118	118	39	74	260
20625-290LW25	123	123	40	77	290
20625-325LW25	128	128	41	80	325
20625-365LW25	133	133	42	83	365
20625-405LW25	138	138	44	86	405
20625-450LW25	143	143	45	89	450
20625-495LW25	148	148	46	92	495
20625-545LW25	153	153	47	95	545
20625-600LW25	158	158	49	98	600

ROUND - EXTRA HIGH PROJECTION					
ORDER N°	DIMENSIONS (mm)				VOLUME (ml)
	A	B	C	D	
20627-195LW25	93	93	47	69	195
20627-225LW25	98	98	49	73	225
20627-260LW25	103	103	51	77	260
20627-300LW25	108	108	53	80	300
20627-345LW25	113	113	55	84	345
20627-390LW25	118	118	57	87	390
20627-440LW25	123	123	59	91	440
20627-495LW25	128	128	61	94	495
20627-555LW25	133	133	63	98	555
20627-615LW25	138	138	65	101	615
20627-685LW25	143	143	67	105	685
20627-755LW25	148	148	69	108	755
20627-835LW25	153	153	71	112	835
20627-920LW25	158	158	73	115	920

ANATOMICAL ROUND AR

textured

PROFILE	BASE	PROJECTION
MODERATE		
HIGH		
EXTRA HIGH		

ANATOMICAL - HIGH PROJECTION					
ORDER N°	DIMENSIONS (mm)				VOLUME (ml)
	A	B	C	D	
20636-190LW25	103	103	43	57	190
20636-220LW25	108	108	44	59	220
20636-250LW25	113	113	46	62	250
20636-285LW25	118	118	48	64	285
20636-320LW25	123	123	49	67	320
20636-360LW25	128	128	51	69	360
20636-405LW25	133	133	53	72	405
20636-450LW25	138	138	54	74	450
20636-500LW25	143	143	56	77	500
20636-550LW25	148	148	58	79	550
20636-610LW25	153	153	59	82	610
20636-670LW25	158	158	61	84	670

ANATOMICAL - MODERATE PROJECTION					
ORDER N°	DIMENSIONS (mm)				VOLUME (ml)
	A	B	C	D	
20635-155LW25	103	103	35	50	155
20635-175LW25	108	108	36	52	175
20635-200LW25	113	113	37	55	200
20635-230LW25	118	118	39	57	230
20635-255LW25	123	123	40	59	255
20635-290LW25	128	128	41	61	290
20635-320LW25	133	133	42	63	320
20635-355LW25	138	138	44	66	355
20635-395LW25	143	143	45	68	395
20635-435LW25	148	148	46	70	435
20635-480LW25	153	153	47	72	480
20635-525LW25	158	158	49	74	525

ANATOMICAL - EXTRA HIGH PROJECTION					
ORDER N°	DIMENSIONS (mm)				VOLUME (ml)
	A	B	C	D	
20637-200LW25	98	98	49	61	200
20637-230LW25	103	103	51	63	230
20637-265LW25	108	108	53	66	265
20637-300LW25	113	113	55	69	300
20637-345LW25	118	118	57	72	345
20637-385LW25	123	123	59	75	385
20637-435LW25	128	128	61	78	435
20637-485LW25	133	133	63	80	485
20637-545LW25	138	138	65	83	545
20637-605LW25	143	143	67	86	605
20637-665LW25	148	148	69	89	665
20637-735LW25	153	153	71	92	735
20637-810LW25	158	158	73	95	810

POLYTECH Health & Aesthetics GmbH
Alzheimer Str. 32 / 64807 Dieburg / Germany
Tel +49 6071 9863 0
Fax +49 6071 9863 30
E-Mail info@polytechhealth.com
www.polytechhealth.com



B-Lite® lightweight breast implants are manufactured in Germany and CE approved.
The products described are covered by one or more US and international patents.
B-Lite® is a registered trademark of G & G Biotechnology Ltd.
Legal manufacturer is G & G Biotechnology Ltd.
© All rights reserved.