



SILICONE SHEETS

Protects the vulnerable skin
around the breast,
during pre-expansion therapy
with EVEBRA-cups.

A PARTNERSHIP BETWEEN:



BAPMEDICAL
MEDICAL CARE FOR THE SKIN



**EVE
PROTECT**

SILICONE SHEETS



Contact us for more information

Aleamed
Duffelsesteenweg 164 a,
2550 Kontich, Belgium.
Tel. +32 3 314 14 01
info@aleamed.be
www.aleamed.eu

EVE-PROTECT SILICONE SHEET

Supportive treatment during pre-expansion therapy in the context of breast augmentation or breast reconstruction using EVEBRA-cups.

- Helps to reduce skin irritation and symptoms of redness and blistering.
- Provides a semi-occlusive (covering) layer on top of the skin, stimulating increased hydration of the skin (epidermis).
- Helps to ensure a good adhesion and sealing around the breast of the EVEBRA-cup.

Product description

- EVE-Protect silicone sheet is a technology for protecting the skin during pre-expansion therapy using EVEBRA-cups (an over-the-breast vacuum device), thereby reducing skin irritation, among other effects.
- This medical, silicone sheet is self-adhesive, washable, and can be used for long periods of time.
- EVE-Protect silicone sheet comes with Cleaning Soap, for cleaning the silicone sheet.

Composition

- Textile layer: polyamide, elastane
- Medical-grade silicone

Characteristics


- EVE-Protect silicone sheet ensures that the skin is protected during pre-expansion therapy using EVEBRA-cups. The sheet helps to reduce skin irritation and symptoms of redness and blistering.
- The silicone sheet also provides a semi-occlusive (covering) layer on top of the skin, stimulating increased hydration of the skin (epidermis).

Indications/applications

EVE-Protect silicone sheet is suitable as supportive treatment during pre-expansion therapy in the context of breast augmentation or breast reconstruction using EVEBRA-cups.

Service life

- Each pack of EVE-Protect contains 1 silicone sheet.



- When used and maintained correctly (according to the instructions for use), the content of each EVE-Protect pack is sufficient for approximately 2 periods of treatment lasting 2 weeks.
- Multiple packs may therefore be necessary for long-term therapy.

Special attention and precautions

- Keep out of reach of children.
- EVE-Protect silicone sheet is suitable for external use only.
- Only use EVE-Protect silicone sheet on healed, closed skin (i.e. after the scab has disappeared) and never on open wounds.
- Never dry EVE-Protect silicone sheet in the sun or near a heat source.
- Temporarily remove EVE-Protect silicone sheet during sports or swimming. Perspiration or contact with water can loosen the EVE-Protect silicone sheet.
- Do not use in case of hypersensitivity to any of the ingredients.
- Do not use if the packaging is damaged on receipt of the product.